



AccuFit® Sentinelle MRI Headrest-Large

REF SHR310



Single use only
Para un solo uso.
À usage unique (jetable) seulement

NON STERILE

LATEX FREE



Medline Industries Inc.
Manufactured exclusively for:
QSUM Biopsy Disposables LLC.
PO Box 2581, Baxter, MN 56425
877-304-4446 Toll Free (218)297-2012 Fax
www.qsum.net Made in Thailand



English- en

1. Open with the elastic ends matching the form of the headrest. Start with the corner, keeping it level pull around to the other end.

Français – French – fr

Ouvrez en faisant les extrémités en élastique correspondre à la forme de l'appui-tête. Commencez par le coin, maintenez en équilibre et tirez le drap jusqu'à en envelopper l'autre extrémité.

Deutsch – German – de

Öffnen Sie das Tuch so, dass die Enden mit Gummizug auf den Enden der Kopfaufgabe zu liegen kommen. Beginnen Sie an einer Ecke, halten Sie das Tuch waagrecht und ziehen Sie es bis zum anderen Ende.

Italiano – Italian – it

Aprire facendo combaciare il bordo dell'elastico con la forma del poggiatesta. Cominciando dall'angolo e mantenendolo stirato, tirare fino all'altro capo.

Español – Spanish – es

Abra por los extremos elásticos haciendo que concuerde la forma del reposacabezas. Comience con la esquina, manténgalo nivelado y tire de él hacia el otro extremo.

Danish –Dansk – dk



QSUM Biopsy Disposables LLC

P.O. Box 2581 Baxter, MN 56425

218.270.3533 Office * 218.297.2012 Fax * 877.304.4446 toll free

* www.qsum.net

* e-mail: Qsumservice@qsum.net

RD12PXS.





Åbn således at elastikenderne matcher hovedstøttens form. Start i hjørnet, hold det lige og træk det til den anden ende.

English- en

2. Pull the drape around the entire headrest.

Français – French – fr

Recouvrez tout l'appui-tête du drap.

Deutsch – German – de

Streifen Sie das Tuch über die gesamte Kopfauflage.

Italiano – Italian – it

Tirare il telo fino a coprire completamente il poggiatesta.

Español – Spanish – es

Tire el paño alrededor de todo el reposacabezas.

Danish –Dansk – dk

Træk afdækningsstykket omkring hele hovedstøtten.



English- en

3. Adjust as necessary.

Français – French – fr

Ajustez selon le besoin.

Deutsch – German – de

Gegebenenfalls anpassen.

Italiano – Italian – it

Aggiustare se necessario.

Español – Spanish – es

Ajuste si es necesario.

Danish –Dansk – dk

Justér om nødvendigt.



**AccuFit and QSUM are registered trademarks of QSUM.*

QSUM Biopsy Disposables LLC

P.O. Box 2581 Baxter, MN 56425

218.270.3533 Office * 218.297.2012 Fax * 877.304.4446 toll free

* www.qsum.net

* e-mail: Qsumservice@qsum.net

RD12PXS.

